



High-Throughput Kit 96x 10 μ mol

Manual

Catalog number: 180484100

(reductively cleavable PEC-Linker RC+)

V1.0 – DRAFT 31.03.2020

Introduction

Belyntic's Peptide Easy Clean (PEC) technology breaks barriers in peptide purification.

Key for a universal application is our novel purification linker, the PEC-Linker.

The target peptides are chemically retained on our activated filter material via robust oxime ligation. Our innovative chemistry allows traceless release of the purified peptides.

This manual provides step-by-step instructions for the PEC purification using our kits, alongside with helpful information and notes to assure a successful use.

Please watch our product video on YouTube (<https://youtu.be/gihpfHr1bfl>) for step-by-step visual guidance. Please also refer to the FAQ on our website [belyntic.com](https://www.belyntic.com).

Now, let's get ready to re-think peptide purification!



Contents

Product Description	5
<i>Kit Content</i>	5
<i>Safety Information</i>	5
<i>Storage Conditions</i>	5
<i>Product Use</i>	5
Preparations	6
<i>Other Chemicals & Consumables</i>	7
<i>Preparation of Buffers</i>	8
Step 1: PEC-Linker Coupling	9
Step 2: TFA Cleavage	10
Step 3: Dissolution & Bead Conditioning	12
Step 4: Immobilization (Catch)	13
Step 5: Washing	13
Step 6: Cleavage & Work-up (Release)	14
Appendix	15
<i>Reaction overview</i>	15
<i>Belyntic's PEC Delta Masses</i>	16
<i>Abbreviations</i>	17
<i>Your Notes</i>	19

Product Description

Kit Content

quantity	label / description	amount
1	PEC filter plate with activated filter material / 50% suspension of Agarose100 provided in 96-well filter plate 400 μ L/well, closed with top- and bottom-sealing-mat	38.4 mL
1	PEC-Linker / RC+	3.2 g
1	Buffer salt / mixture of citric acid-sodium carbonate	5.6 g
1	Blocking agent / L-cysteine	1.0 g
1	Reducing agent / Dithiothreitol	5.0 g
1	TFA collection plate / deep-well plate to collect crude peptides, perform ether precipitation and centrifugation closed with top-sealing-mat	
1	Peptide collection plate / perforated deep-well plate to collect purified peptides and perform freeze-drying in the plate	

Safety Information

Please read the Safety Data Sheets, available at [belyntic.com](https://www.belyntic.com) before use.

Storage Conditions

Upon receipt, store the kit components in the fridge at 2-8°C. The agarose beads should not be frozen or dried under reduced pressure. Buffers (cf. page 8) can be stored for 2 weeks in the fridge at 2-8°C.

Product Use

For Research Use Only. The use in or for clinical or diagnostic products, applications, or services is excluded.

Preparations

Purification Considerations

- PEC-Linker RC+ will be installed at the *N*-terminus of the peptide
- reductive conditions will be applied in the purification process
- azide and disulfide groups in the peptide might be reduced

Precautions During Peptide Synthesis

- **acetylation after each amino acid coupling (capping) step is key to ensure a successful PEC purification**

efficient conditions for blocking unreacted amino groups are 2 M acetic anhydride (20 vol%) with either 2 M pyridine (18 vol%) or 2 M 2,6-lutidine (21 vol%) in DMF for 5 min¹

- **special attention must be paid to possible side reactions during SPPS**

- minimize by-products within the fulllength chain (e.g. aspartimide)
- minimize by-products during TFA cleavage¹, ideally use Fmoc-Lys(Trt)-OH to prevent *t*Bu on Lysine² and Fmoc-Cys(StBu)-OH to prevent alkylation on -SH during TFA cleavage
- use Ramage- or Sieber-amide resins for C-terminal peptide amides and chlorotrityl resin for C-terminal peptide acids (avoid Wang-linker)
- N-terminal glutamine bearing sequences may show 10-20% of pyroglutamate after purification, contact us for support if you want to purify N-terminal glutamine peptides

- **avoid aldehydes and/or ketones in ether and solvents**

the free oxy-amine on PEC-linker will readily react with traces and contaminations of ketones and/or aldehydes that will deactivate it

¹ O. Reimann, O. Seitz, D. Sarma, R. Zitterbart, *J. Pept. Sci.* **2019**, 25, e3136.

² J. Pawlas, U. Törnvall, *Proceedings of the 33rd European Peptide Symposium* **2014**, 108-109.

Other Chemicals & Consumables

Reagents

- ethyl-2-cyano-2-(hydroxyimino) acetate (Oxyma) or 1-hydroxybenzotriazol (HOBt) as coupling reagent
- scavengers for TFA cleavage cocktail, e.g. TIS, EDT, thioanisol, phenol
- guanidinium chloride (GdmCl) - in case of technical grade, filter the solution prior to use
- sodium chloride (NaCl)

Acids / Bases

- *N,N*-diisopropylethylamine (DIEA)
- trifluoroacetic acid (TFA)
- sodium hydrogen carbonate (NaHCO₃)

Solvents

- *N,N*-dimethylformamide (DMF), ≥ 99.5% (peptide synthesis grade)
- dimethyl sulfoxide (DMSO), ≥ 99.5% (synthesis grade)
- dichloromethylene (CH₂Cl₂), ≥ 99.5% (synthesis grade)
- acetonitrile (MeCN), ≥ 99.5% (synthesis grade)
- deionized water
- ethanol (EtOH), denatured
- suitable ether for peptide precipitation (Et₂O and n-heptane 3:1), ≥ 99.5% (synthesis grade, acetone-free)
- 1,1,1,3,3,3-hexafluoro-propan-2-ol (HFIP), if using synthesis grade (≤99%) please reflux 1 h with 1 w% L-cysteine to remove trifluoroacetone, distill then to get pure HFIP



Preparation of Buffers

□ **Buffer#1**

dissolve 5.6 g of the supplied buffer salt slowly in 200 mL water (final pH = 4.5, caution: strong CO₂ evolution); To prepare more Buffer#1 see below³

□ **Buffer#2 (7 M GdmCl in Buffer#1)**

dissolve 13.4 g of GdmCl in 20 mL Buffer#1

□ **Wash A: DMSO with 0.9 M GdmCl**

dissolve 13 g of GdmCl in 150 mL of DMSO

□ **Wash B: EtOH/H₂O (7:3) with 0.1 M NaCl**

dissolve 1168 mg of NaCl in 60 mL water **first**, then add 140 mL EtOH

□ **NaHCO₃ pH 8:**

dissolve 2.5 g of NaHCO₃ in 50 mL water, to give a pH of 8

Equipment & Processing

please don't hesitate to contact us, if in doubt about suitable equipment

- shaking device for deep-well plates, ideally with 800 rpm
- centrifuge with rotor for deep-well plates
- Belyntic's Acrylic Vacuum Manifold, catalog number: 180484150. Other manifolds cannot be recommended since they may produce cross-contaminations
- 30-300 µL multichannel pipette (8- or 12-channel) or dispenser with 8- or 12-channel manifold exit

³Buffer#1: dissolve 19.2 g Citric Acid (anhydrous, or 21.0 g of monohydrate) together with 8.75 g Na₂CO₃ (anhydrous) slowly in 1 L water.

Step 1: PEC-Linker Coupling

Note: The following procedure describes manual coupling of PEC-Linker RC+ after peptide assembly in a synthetic filter plate. The PEC-Linker can also be coupled to the peptides on a synthesizer; higher volumes might be necessary due to synthesizer dead volume. The base (DIEA or *N*-methylmorpholine) must be added to the pre-mixed solution of the PEC-Linker with Oxyma (or HOBt) before coupling.

- swell synthetic resin in DMF for 15 min and remove DMF before PEC-Linker addition
- weigh in the amount of PEC-Linker and Oxyma or HOBt (Table 1) and combine in a suitable vessel

Table 1: PEC-Linker coupling for 96x 10 μ mol peptide

PEC-Linker RC+ (>4 eq.)	3200 mg	4255 μ mol
Oxyma/ HOBt (>6 eq.)	907/862 mg	6383 μ mol
add DMF to have ~13.8 mL	10630 μ L	0.3 M PEC-Linker
DIEA (>6 eq.)	1116 μ L	6383 μ mol

- add DMF to dissolve PEC-Linker and Oxyma or HOBt (Table 1)
- then add DIEA (Table 1, on a synthesizer DIEA or *N*-methylmorpholine should be added directly in a pre-activation vessel before coupling)
- quickly add 150 μ L of the yellow solution to each well

Note: Other common solvents can be used (e.g. CH_2Cl_2 , DMSO or NMP)).

- couple for ≥ 3 h (shaking is not necessary; overnight is possible)
- wash synthesis resin with 3x ≥ 200 μ L/well DMF and 3x ≥ 200 μ L/well CH_2Cl_2

Note: You may pause here and store the linker-modified peptide on synthesis resin below 8°C.

Step 2: TFA Cleavage

- choose and prepare proper cleavage cocktail for your peptides (500 μ L/well = 48 mL)

By-products caused by re-addition of protecting groups cannot be removed with PEC, please use proper scavengers: Reagent B1 TFA/H₂O/PhSMe/EDT/TIS (83:5:5:5:2) or a non-odorous cocktail of TFA/H₂O/DTT/TIS (88:4:6:2) are recommended.

- remove the sealing-mat, put the TFA collection plate in a vacuum block from your synthesizer and place the synthesis plate on the manifold
- in a fume hood, add TFA-cleavage cocktail to each well of the synthesis plate and carry out cleavage for 2 h (no shaking); use reduced pressure at the end to collect cleavage solution
- take the TFA collection plate that is containing the peptides in TFA cocktail out of the manifold and precipitate the peptides by addition of 1.2 mL **cold** (≤ -20 °C) Et₂O/n-Heptane (3:1) to each well

Aldehyde- and ketone-free ether is crucial during this step, in doubt use THF/n-Heptane (1:1) as an alternative.

- close the plate with the sealing-mat (**make sure that it is closed properly**), shake for 30 seconds and put it in the freezer (≤ -20 °C) for ≥ 10 min
- centrifuge the plate for ≥ 4 min at ≥ 4000 rpm and discard the ether supernatants by carefully turning the plate upside down. The precipitates will stay in the TFA collection plate.
- put the TFA collection plate back in the vacuum manifold
- place the synthesis plate on the manifold and elute each well with 2x 200 μ L TFA/H₂O (95:5), use reduced pressure to decrease elution time
- take the TFA collection plate out of the manifold

If you have methionine containing peptides, please add 7 μL TMSBr and 8 μL EDT to each well containing such peptides and shake the plate gently (300 rpm) for 15 min.

- precipitate the peptides again by addition of 1.2 mL **cold** ($\leq -20\text{ }^\circ\text{C}$) $\text{Et}_2\text{O}/n\text{-Heptane}$ (3:1) to each well
- close the plate with the sealing-mat (**make sure that it is closed properly**), shake for 30 seconds and put it in the freezer ($\leq -20\text{ }^\circ\text{C}$) for ≥ 10 min
- centrifuge the plate and discard the ether supernatant as described above
- let the ether precipitates dry by putting the plate uncovered in a fume hood for ≥ 30 min
- cover the TFA collection plate with a paper tissue fixed with rubber bands and dry the precipitates carefully in an exicator with a membrane pump (**after** drying in a fume hood)
- clean the sealing-mat with water

Note: You may pause here and store the crude linker-modified peptides as a dried solid below $8\text{ }^\circ\text{C}$.

If crude analysis is performed, check delta mass list in Table 2 (cf. page 16). On some analytical LC systems, the linker-modified peptide cannot be detected or shows lower purities. Dilute your sample with a 50 mg/mL solution of methoxyamine hydrochloride (CAS: 593-56-6), that will prevent analytical artifacts due to aldehydes or ketones. However, if you still don't see linker-modified peptide, please continue the purification process described in this manual as the desired peptide will be visible again after PEC-Linker cleavage.

Step 3: Dissolution & Bead Conditioning

- first add 450 μL /well DMSO and close the plate with the sealing-mat

Note: For very hydrophobic peptides HFIP can be used (make sure it is free of trifluoroacetone, solvent section page 7). Other solvents, such as DMF or MeCN up to 90 vol% in Buffer#1, may be used. Ketone- or aldehyde-containing solvents are not tolerated.

- shake heavily for ≥ 1 h (overnight possible)

Please make sure to use a shaker that is suitable for deep-well plates and use a shaking speed of 800 rpm.

Before removing sealing-mats, centrifuge the plates for 1 min at ≥ 4000 rpm to prevent cross-contaminations

If you have methionine containing peptides and you want to use DMSO, please use DMSO purged with inert gas (N_2 or Ar) and try not to exceed dissolution time of 1 h and immobilization time of 4 h.

- if all peptides are dissolved (inspection with flashlight), please add 50 μL /well Buffer#2; **don't add Buffer#2 when using HFIP as solvent**
- take the PEC filter plate that is prefilled with agarose suspension provided in your High-Throughput Kit, remove top- and bottom-sealing-mat, place it in the vacuum manifold and remove the supernatants
- wash 3x 300 μL /well with water and 3x 300 μL /well with Buffer#1 (2.4 mL/8-wells, 3.6 mL/12-wells, 90 mL in total each)
- wash once with 300 μL /well DMSO **or again with 300 μL /well water when using HFIP as peptide solvent**
- close the PEC filter plate with the bottom-sealing-mat and add 200 μL /well of the solvent used for peptide dissolution

Step 4: Immobilization (Catch)

- add the dissolved peptides to their respective wells
- apply the top-sealing-mat and immobilize peptides by shaking for ≥ 4 h (overnight possible, when not dealing with methionine containing sequences)
- remove the sealing-mats from the PEC filter plate and drain the immobilization supernatants by usage of the vacuum manifold
- wash 3x 300 μ L/well with DMSO (2.4 mL/8-wells, 3.6 mL/12-wells)

Note: You may monitor immobilization efficiency by LC or MS analysis of the supernatant of single wells prior to addition of Blocking Buffer.

- prepare Blocking Buffer (2 w% solution of L-cysteine): dissolve all L-cysteine (1 g) provided in your High-Throughput Kit in 50 mL Buffer#1
- close the PEC filter plate with the bottom-sealing-mat and add 500 μ L/well of Blocking Buffer
- apply top-sealing-mat and shake for 30 min to block unreacted aldehyde groups and to liberate imines
- remove the sealing-mats, place the PEC filter plate on the vacuum manifold and remove the Blocking Buffer

Step 5: Washing

Each wash is 3x 300 μ L/well (2.4 mL/8-wells, 3.6 mL/12-wells, 90 ml in total)

- **wash A:** DMSO with 0.9 M GdmCl
- **wash B:** EtOH/H₂O (7:3) with 0.1 M NaCl

Note: You may pause here by storing the immobilized peptides after one additional H₂O wash in the fridge at 2-8 °C.

Step 6: Cleavage & Work-up (Release)


- prepare Reducing Agent: dissolve 5.0 g DTT in 50 mL 5 w% aq. NaHCO₃ (0.6 M; pH 8)
- close the PEC filter plate with the bottom-sealing-mat and add 500 µL/well MeCN
- add 500 µL/well of the Reducing Agent, apply top-sealing-mat and shake for 30 min, shake for 90 min when Cys(StBu) has been used
- wash 3x 300 µL/well with water, 3x 300 µL/well with MeCN and discard to remove excess DTT

From this point on, please collect all filtrates in the perforated peptide collection plate, since they hold your product.

- put the peptide collection plate in the vacuum manifold and place the PEC filter plate on the manifold
- add 200 µL/well H₂O/MeCN/TFA (90:8:2)
- let the plate stand in the vacuum manifold for ≥ 4 h (overnight possible)
- collect the peptides by reduced pressure and further elute peptides:
1st: 600 µL/well H₂O/MeCN (7:3) + 0.1% TFA
2nd: 500 µL/well H₂O/MeCN (3:7) + 0.1% TFA
(use reduced pressure to decrease elution times)
- take the peptide collection plate out of the vacuum manifold and take your analytical samples
- cover the plate with a paper tissue fixed with rubber bands, freeze in liquid nitrogen **for ≥ 30 min** and lyophilize the peptides

Appendix

Reaction overview





Belyntic's PEC Delta Masses

Table 2: Relevant delta masses during the purification process.

+412.03/414.03 (MH ²⁺ : +206.52, MH ³⁺ : +137.67)	PEC-Linker RC+ on peptide
+ 40.05 (MH ²⁺ : +20.03, MH ³⁺ : +13.33)	acetone oxime of peptide-linker conjugate, +453.08 Da with PEC-Linker RC+
+ 26.03 (MH ²⁺ : +13.02, MH ³⁺ : +8.68)	acetaldehyde oxime of peptide-linker conjugate, +439.06 Da with PEC-Linker RC+
+ 15.99 (MH ²⁺ : +8.00, MH ³⁺ : +5.33)	oxidation of methionine residues
+ 12.01 (MH ²⁺ : +6.01, MH ³⁺ : +4.00)	formaldehyde oxime of peptide-linker conjugate, +425.04 Da with PEC-Linker RC+
+56.07 (MH ²⁺ : +28.04, MH ³⁺ : +18.69)	reattached tBu; thiols (e.g. EDT) should be used in the TFA cleavage cocktail (make sure they are not oxidized)
+106.05 (MH ²⁺ : +53.03, MH ³⁺ : +35.35)	para-hydroxybenzyl cation, generated during TFA cleavage; use RAM-linker or add EDT and 1,3-Dimethoxybenzene during cleavage
+252.06 (MH ²⁺ : +126.03, MH ³⁺ : +84.02)	non-cleaved or reattached Pbf; treat purified peptide with Reagent K or another thioanisole-containing TFA cocktail for 3-4 h.

Abbreviations

AcOH	acetic acid
Buffer#1	immobilization buffer
Buffer#2	7 M guanidinium chloride in immobilization buffer
CAS	Chemical Abstracts Service
DIEA	<i>N,N</i> -Diisopropylethylamine
DMF	<i>N,N</i> -Dimethylformamide
DTE	dithioerythritol (Cleland's reagent)
DTT	dithiothreitol (Cleland's reagent)
EDT	ethanedithiol
Et ₂ O	diethyl ether
EtOH	ethanol
GdmCl	guanidinium chloride
HFIP	1,1,1,3,3,3-hexafluoro-propan-2-ol
HOBt	1-hydroxybenzotriazole
LC	liquid chromatography
MeCN	acetonitrile
MS	mass spectrometry
MTBE	<i>t</i> -butyl methyl ether
NMP	<i>N</i> -methyl pyrrolidone
Oxyma	ethyl-2-cyano-2-(hydroxyimino) acetate
Pbf	pentamethyldihydrobenzofuran-5-sulfonyl
PEC	Peptide Easy Clean
PhSMe	Thioanisol
RAM	modified Rink-Amide
SDS	safety data sheet
SPPS	solid phase peptide synthesis
<i>t</i> Bu	<i>t</i> -butyl
TFA	trifluoro acetic acid
THF	tetrahydrofuran
TIS	triisopropylsilane
TMSBr	trimethylsilyl bromide



Technical Support

Please watch our product video on YouTube for visual step-by-step guidance (<https://youtu.be/gihpfHr1bfl>). Please also refer to the FAQ on our website <http://belyntic.com>.

For further technical support, submit your questions directly via mail to support@belyntic.com or call +49 30 8104-1113.

Safety Data Sheets (SDS)

Safety Data Sheets (SDS) are available on our website <http://belyntic.com>.

Purchaser Notification

Please read the Terms of Sale available on our website <http://belyntic.com>.

Contact

Belyntic GmbH

Richard-Willstaetter-Str. 11

12489 Berlin, Germany

Phone: +49 30 8104-1113

Email: support@belyntic.com

belyntic.com